



A new approach to mask selection

Distinct categories for every patient preference



An easier way to help patients choose the right mask

Helping your patients choose the right mask at setup is one of the most important parts of CPAP therapy success. ResMed mask categories make it easy to understand key benefits, and can serve as a conversation starter between you and your patients. Every ResMed mask is designed to deliver effective therapy and a comfortable night's sleep, but each category offers something a little different.

Minimalist category

Small. Light. Simple.

Freedom category

Sleep in any position.

Ultra Soft category

Patented memory foam offers unique comfort.

Universal Fit category

Classic designs that fit nearly every face.



For many patients, their primary driver for mask choice is simple: less mask on their face. This category features masks that are small, lightweight and easy to use. Minimal facial coverage is also appealing to anyone who's claustrophobic.

This category features masks with a top-of-head connection, which keeps the system's tubing out of sight and out of the way. These masks help active sleepers comfortably move around at night and sleep in whichever position they choose, including on their stomach.

Mask discomfort is the number one reason CPAP users fall off therapy¹ and the Ultra Soft category is designed to change that. With its patented memory foam cushion, AirTouch™ F20 delivers a unique level of comfort and softness.

With a fit range of 97% for AirFit™ F20²⁻⁴ and 99% for AirFit N20⁵, you can always reach for these masks with confidence. They're also designed to perform across a wide range of therapy pressures.



AirFit™ N30



AirFit™ F30



AirFit™ P10



AirFit™ F30i



AirFit™ N30i



AirFit™ P30i



AirTouch™ F20



AirFit™ F20



AirFit™ N20



Contact your ResMed sales rep
to learn more.



1 ResMed AirTouch™ F20 rescue study: 30-day assessment of 30 participants who intended to quit CPAP therapy with prescribed mask at baseline. Jan. 15, 2018 to Feb. 22, 2018 in San Diego, CA.
2 ResMed AirFit F20 internal Australian fitting study of existing CPAP patients conducted between March–April 2016.
3 ResMed AirFit F20 internal USA fitting study of existing CPAP patients conducted April 2016.
4 ResMed AirFit F20 internal EU and APAC fitting study of existing CPAP patients conducted June 2016.
5 ResMed AirFit N20 internal international fitting study of new and existing CPAP patients, conducted Nov 2015.
Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA. (858) 836-5000 or (800) 424-0737 (toll free). See ResMed.com for other ResMed locations worldwide. AirFit and AirTouch are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. For patent and other intellectual property information, see ResMed.com/ip.
© 2019 ResMed. 10111432/3 2019-11

ResMed.com